

## GE Healthcare Imagescanner III



The GE Image Scanner III is used for acquiring images from two-dimensional (2-D) electrophoresis gels and blots using the included LabScan software

Rating: Not Rated Yet

**Price**

Sales price \$3,100.00

[Ask a question about this product](#)

Manufacturer [GE Healthcare](#)

### Description

**!!!! - RECENTLY SOLD - !!!!**

The GE Image **Scanner** III is a modified version of the Epson Expression 10000XL **Scanner**, it is used for acquiring images from two-dimensional (2-D) electrophoresis gels and blots using the included LabScan software, comes with a computer with LabScan and Epson Scan software installed, monitor, keyboard and mouse, user manuals included, 100-240V, 50-60Hz, 1.2-0.7A.