


Gilson Minipuls 3 Peristaltic Pump with 4-Channel Head



Flow Rate: 0.3 μ L/min to 30 mL/min

Rating: Not Rated Yet

Price

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Gilson](#)

Description

The Gilson MINIPULS® 3 Peristaltic Pump with 4-Channel Head is a compact, high-performance pump designed for precise, low-pulse liquid handling in laboratory settings, delivering a variety of biological and chemical fluids with minimal shearing using interchangeable pump heads and microprocessor-controlled stepper motor technology.

Key Specifications:

Flow Rate: 0.3 μ L/min to 30 mL/min (standard 4-channel head); up to 220 mL/min with high-flow 4-channel head

Channels: 4 channels, supporting simultaneous delivery through multiple tubes.

Pressure: Up to 0.5 MPa (standard flow) or 0.3 MPa (high flow)

Control: Microprocessor-controlled with adjustable speed stepper motor, keypad, and LCD for start/stop, flow direction, rate adjustment, and priming

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com