

Heidolph Basis Hei-VAP HL with Vertical Coil Condenser



Temperature range (°C): 20 - 210; 4.5 L Capacity.

Rating: Not Rated Yet

Price

Sales price \$3,500.00

[Ask a question about this product](#)

Manufacturer [Heidolph](#)

Description

Rotary Evaporator with Vertical Coil Condenser.

Model: Heidolph Hei-VAP Value Digital

Specifications:

Temperature range (°C): 20 - 210

Rotation speed: 10 - 280 1/min¹

Heating capacity: 1400 W

Cooling surface: 1400 cm²

Bath capacity: 4.5 L

Lift distance: 155 mm

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com