

Heidolph Multi Reax Rotator



Speed Range: 150–2,500 rpm

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Heidolph](#)

Description

The Heidolph Multi Reax is a versatile vibrating test tube shaker designed for efficient mixing of samples in tubes and vials of various diameters, ideal for scientific research, chemistry, microbiology, and pharmacology applications.

Key Specifications:

Speed Range: 150–2,500 rpm, continuously adjustable via digital display.

Vibration Orbit: 3 mm, suitable for rigid to smooth dispersion of viscous media and solids.

Tube Capacity: Includes carousel attachments for up to 26 tubes (10–16 mm diameter, max 160 mm height) and 12 tubes (16–32 mm diameter, max 120 mm height).

Timer: Digital process timer for 1–999 minutes with acoustic alarm and automatic shut-off.

Display: Four-digit LED showing elapsed time, rotational speed, and start/stop indicators.

Shakers/Stirrers/Hotplates: Heidolph Multi Reax Rotator

Operation Modes: Continuous shaking or timer-controlled with automatic shut-off.
Motor: Sparkless, maintenance-free with self-reset function for overheating protection.
Power Supply: 230/240V or 115V, 50/60Hz

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com