


Hettich Rotina 35



Hettich Rotina 35 high-performance benchtop centrifuge w/2 rotors

Rating: Not Rated Yet

Price

Sales price \$3,250.00

[Ask a question about this product](#)

Manufacturer [Hettich](#)

Description

!!!SOLD!!!

This Hettich Rotina 35 high-performance benchtop centrifuge comes w/2 rotors, a 10 by 50 ml swinging bucket one (to 5,000 rpm) and a two-place titerplate rotor (to 4,000 rpm). 1. Automatic lid locking, auto lid lock. 2. Frequency-controlled induction drive. 3. Complete running programs can be stored and recalled. Can reach up to 15,000 rpm without fixed angle rotor. Very, very clean. Manufactured in 2002.

Please call us at (650) 871-5707 to check availability of this product!