



## Autoclaves: Hirayama HG-80 Upright Autoclave

---

Programmable Auto-Start for running a cycle from 1 to 99 hours later.

Large (15 liters)/ Efficient Steam Condenser ensures no steam exhausts into the lab.

Over-Flow Lamp indicates that the drain bottle is full and ready to be emptied. Built-in Cooling Fan for faster post-sterilization cooling and shorter completion time.

Selectable Rate of steam venting after sterilization.

Programmable Time for air removal before sterilization.

Twelve Programs Cycles can be stored.

Large Control and Display Panel improves functionality and shows the progress of the sterilization cycle.

Four Program Cycles And Three SubPrograms 1) Liquid cycle, 2) Agar (liquid with warming) cycle, 3) Solid cycle; and 4) Agar dissolution and warming. Each of these four programs has 3 sub-programs in which sterilization parameters can be stored, providing a total of 12 stored program cycles.

Other Specifications:

Temperature Range: 105C to 135C

Sterilization Timer : 1 to 300 minutes

Media Warming/ Holding Temperature: Programmable at 45C to 60C to keep media warm up to 24 hours Agar Dissolution & Warming Protocol Initially at 60C-100C for 1-60 minutes, then at 45C to 60C for up to 24 hrs.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)