

## Hygiena SystemSURE Plus Meter



Measurement Method: ATP bioluminescence with solid-state photodiode detection

Rating: Not Rated Yet

**Price**

Sales price \$1,150.00

[Ask a question about this product](#)

### Description

The Hygiena SystemSURE Plus is a compact, handheld ATP luminometer used for rapid hygiene verification in food safety, healthcare, water treatment, and sanitation programs. Used with Hygiena ATP swabs, the unit provides quick quantitative pass/fail results in seconds, helping verify surface cleanliness and validate sanitation procedures in the field or lab.

### Specifications

Measurement Method: ATP bioluminescence with solid-state photodiode detection

Display: Backlit LCD screen with simple menu-driven navigation

Memory Storage: Up to 250 user locations / test plans (model dependent)

Calibration: Internal calibration – no external calibration required

Results Time: < 15 seconds per test

Data Interface: USB connection for data export and software integration

Power: AA alkaline or rechargeable batteries (typically >8 hours of use per set)

Compatible Swabs: UltraSnap, SuperSnap, AquaSnap, and other Hygiena ATP tests

Size & Weight: Small, lightweight handheld unit designed for field portability

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-

quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax

[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)