

IEC Centra CL3 Centrifuge



IEC Centra CL3 Centrifuge (Rotor not included)

Rating: Not Rated Yet

Price

Sales price \$1,900.00

[Ask a question about this product](#)

Manufacturer [IEC](#)

Description

The IEC Centra CL3 Centrifuge:

- Store and lock up to 99 protocols
- Speed and accuracies to the nearest 50 rpm as well as automatic speed calculation of the g-force
- Temperature with 1 °C increments
- Capable of reaching 4 °C from 20 °C in ~7 minutes with rapid condition means you always spin at the required temperature.

The CL3 series accommodates up to 1 liter of sample (4 x 250mL) accepting a wide variety of popular tubes including 15 mL, 50 mL, 225 mL conical, Vacutainer® tubes to 15 mL (16 x 125mm), microtubes, microplates (standard and deepwell), microscope slides, and an assortment of filtration devices.

New Quick Release Rotor

The Thermo IEC 243 Swinging Bucket Rotor, utilizes convenient sealed Aerocarrier Adapters which accommodate sample volumes ranging from 0.25 mL to 250 mL.

****Centrifuge does not include rotor****