

IEC Micromax RF



IEC Micromax RF refrigerated digital microprocessor controlled microfuge

Rating: Not Rated Yet

Price

Sales price \$1,650.00

[Ask a question about this product](#)

Manufacturer [IEC](#)

Description

SOLD

This IEC Micromax RF refrigerated digital microprocessor controlled microfuge has a 24 by 1.5-2.0 ml rotor, can get to speeds of 15,000 rpm, a temperature range of -9 to 40C and has a 1-99 minute timer (smaller increments for shorter runs). Has a brushless motor and is very quiet. Very light (70 lbs) for a refrigerated unit. Manufactured in 2000.

Please call us at (650) 871-5707 to check availability of this product!