


## IKA Plate S001 Digital Heating Magnetic Stirrer





Temperature Range: Up to 310°C; Speed Range: Up to 1500 RPM

Rating: Not Rated Yet

**Price**

Sales price \$425.00

[Ask a question about this product](#)

Manufacturer [IKA](#)

Description

IKA PLATE S001 Digital Heating Magnetic Stirrer

Rapid and uniform heating along with precise digital temperature control. For a wide range of applications including titrations, sample preparations and chemical syntheses.

Specifications:

Platform: 5" Diameter, Aluminum Alloy

Temperature Range: Up to 310°C

Speed Range: Up to 1500 RPM

Weight Capacity: 55 lb

Stirring Volume: 20 L

Display: LCD

Ports: PT1000, ETS-D5, and ETS-D6

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)