


IKA Vibrax VXR Basic



wide speed range (0 - 2200 rpm)

Rating: Not Rated Yet

Price

Sales price \$550.00

[Ask a question about this product](#)

Manufacturer [IKA](#)

Description

Optoelectronically controlled small shaker with a very wide speed range (0 - 2200 rpm).

Features:

- Suitable for continuous operation
- New design and improved drive system
- Circular shaking motions
- Slow speeds are well maintained

Attachments are interchangeable

Specifications:

Type of movement: orbital
Shaking stroke: 4 mm
Permissible shaking weight (incl. attachment): 2 kg
Motor rating input: 35 W
Motor rating output: 13.5 W
Permissible ON time: 100 %
Speed min (adjustable): 100 rpm
Speed range: 0 - 2200 rpm
Speed display: scale
Operating mode: continuous operation
Dimensions (W x H x D): 157 x 130 x 247 mm
Weight: 6.1 kg
Permissible ambient temperature: 5 - 40 °C
Permissible relative humidity: 80 %
Voltage: 100 - 120 V
Frequency: 50/60 Hz
Power input: 45 W

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com