

## Integra Biosciences MINI 96 Automated Pipettor 12.5µL



Manufactured in 2024.

Rating: Not Rated Yet

**Price**

Sales price \$7,750.00

[Ask a question about this product](#)

### Description

Model no. 4801 The MINI 96 is a compact and portable semi-automated pipettor capable of pipetting into 96- and 384- shallow and deep well microplate formats. Unit features the 12.5µL pipetting head for volumes between 0.5-12.5µL.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)