

Kuhner ISF-4-W Large Capacity Incubator Shaker



The Climo Shaker (model ISF-4-V) is a floor model incubator/shaker featuring dual side by side doors with double pane glass.

Rating: Not Rated Yet

Price

Sales price \$19,500.00

[Ask a question about this product](#)

Manufacturer [Kuhner](#)

Description

!!!SOLD!!!

The unit is equipped with four stainless steel shelves each with variable speed orbital shakers able to range from 30/300 RPM. Digital displays on the outer panel control the shaking feature as well as the temperature inside the unit, which ranged from 15 to 60 degrees Celsius. The Climo Shaker operates on 200V single phase power.