

Lab-Line 4625 Titer Plate Shaker



Lab-Line 4625 Titer Plate Shaker holds up to 4 plates and gets up to 1100 rpm

Rating: Not Rated Yet

Price

Sales price \$450.00

[Ask a question about this product](#)

Manufacturer [Lab-Line](#)

Description

This Lab-Line 4625 Titer Plate Shaker holds up to 4 plates and gets up to 1100 rpm. Electrical timer can be set for continuous operation or timed cycles up to 5 minutes. White, non-slip rubber mat increases sample visibility. Multiple units available.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com