

## LabArmor Bead Bath



With 12L of beads.

Rating: Not Rated Yet

**Price**

Sales price \$995.00

[Ask a question about this product](#)

**Description**

Model 74300-714

Produces a thermal envelope that is highly uniform across the entire bath.

With 12L of beads.

**Specs:**

Temperature Uniformity: +/- 1.0°C;

Temperature Range: 5°C above ambient to 80°C

Exterior Dimensions WxDxH in. (cm): 15.9 x 18.5 x 8.5 (40.5 x 40.5 x 20.9); Interior Dimensions WxDxH in. (cm): 12 x 12 x 6 (30.5 x 30.5 x 15)

Electrical requirements: 4.5A, 100-120V, 50/60Hz

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)