

Labnet Prism R Centrifuge



Up to 13,500 rpm or 17,135 × g

Rating: Not Rated Yet

Price

Sales price \$2,650.00

[Ask a question about this product](#)

Manufacturer [Labnet](#)

Description

A refrigerated microcentrifuge that combines a compact footprint with a powerful brushless motor and rapid cooling capability, reaching 4°C in under 8 minutes. It features an intuitive LCD control panel for precise settings of speed (up to 13,500 rpm or 17,135 × g), time, and temperature, along with safety functions like imbalance detection and vibration-free operation. The 24-place aluminum rotor supports 1.5/2.0 mL tubes directly, accommodates smaller tubes with adapters, and allows for an optional attachment for 0.2 mL PCR tubes and strips.

Specifications

Max RCF: 17,100

Centrifuges: Labnet Prism R Centrifuge

Spins these tubes: 1.5-2.0ml, .2mL strip tubes

Maximum rotor capacity: 0.2 PCR strip tube x 2, 24 x 1.5/2.0mL

Timer: 0.5 to 99 minutes or continuous, "Quick" button for momentary operation

Voltage: 120V

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com