

## Leica CM1950 Cryostat



Maximum specimen size: 50 x 80 mm

Rating: Not Rated Yet

**Price**

Sales price \$22,500.00

[Ask a question about this product](#)

Manufacturer [Leica](#)

### Description

Model #: CM1950

The cryostat for laboratories with a large number of specimens and varying specimen types. Provides sufficient space to organize your tools, freezing media and glass slides.

- Width (with handwheel): 835 mm
- Depth (cabinet only): 850 mm
- Height (total): 1215 mm
- Working height (armrest): 1025 mm

- Weight (w/ motor and suct.): 193 kg
- Operating temperature range: 18°C to 35°C
- Maximum specimen size: 50 x 80 mm
- Electric coarse feed:
  - Slow: 300 ?m/s, in 20 ?m steps
  - Fast: 900 ?m/s

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)