

Lonza FlashGel Device Pack



Lonza 57068 FlashGel with Lonza 57025 Flash Gel

Rating: Not Rated Yet

Price

Sales price \$975.00

[Ask a question about this product](#)

Description

The FlashGel™ Device Pack provides the key components for researchers to quickly begin using and imaging FlashGel™ DNA Cassettes. FlashGel™ DNA Cassettes are the ideal sample screening tools to check up to 32 PCR or restriction fragments quickly without having to plan your day around agarose gels. Monitor and document gel runs right at the bench, without UV light. The highly sensitive system allows a 5-20X reduction in DNA levels--saving both time and money. FlashGel™ Cassettes contain precast, prestained agarose gels and buffer--no need for gel preparation, buffer addition or gel staining.

Specifications:

Standard: Metric

Manufacturer Name: Lonza

Serial Number: 170321044 112187

Model: 57068

Voltage: 100-240 V

Frequency: 50/60 Hz

Recommended Packaging Form: 10C (10X10X10)

Dimensions Overall: 15 x 6 x 4 IN - 2 lbs / 381.00 x 152.40 x 101.60 mm - 0.91 kg

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com