

Luminex 100



Luminex 100 (Bio-Rad Bio-Plex) system w/main unit, XY platform, and SDS module. With copy of MiraiBio MasterPLex software

Rating: Not Rated Yet

Price

Sales price \$16,500.00

[Ask a question about this product](#)

Manufacturer [Luminex](#)

Description

!!!! - RECENTLY SOLD - !!!!

This Luminex 100 (Bio-Rad Bio-Plex) system w/main unit, XY platform, SDS module, and copy of MiraiBio MasterPLex software is guaranteed installable by a Luminex engineer.

The Luminex 100 System can deliver information accurately, inexpensively and in real-time enhancing your operational efficiency. This compact lab analysis system includes the Luminex 100, the XY Platform, and the SD. It utilizes xMAP® technology which enables you to simultaneously assay up to 100 analytes in a single well of a microtiter plate, using very small sample volumes. The system delivers fast and cost-effective bioassay results on many assay formats including nucleic acid assays, receptor-ligand assays, immunoassays and enzymatic assays.