

## Manostat Simon Varistaltic Pump



This unit is perfect for moving liquids and gasses, transferring chemicals out of carboys, and moving gas through flow systems.

Rating: Not Rated Yet

**Price**

Sales price \$375.00

[Ask a question about this product](#)

Manufacturer [Manostat](#)

### Description

This unit is perfect for moving liquids and gasses, transferring chemicals out of carboys, and moving gas through flow systems. It can also be used to generate pressures up to 20 psi. Variable speed control capable of flow rates from 1 ml/minute to 3450 ml/minute. Precise electronic controls with analog allow infinite flow rate selection adjustment. Reversible speed direction for purging of transfer lines. Self-priming pump head accepts any flexible tubing with 1/16" wall and I.D. 1/32" to 5/16". Dimensions: 26L x 23D x 13H cm. Power: 115 volt, 50/60 Hz.

*Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.*

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)