

Vacuubrand MD1C Vacuum Pump



.88 cfm/ 1.5Torr, (2mbar)

Rating: Not Rated Yet

Price

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [Vacuubrand](#)

Description

VACUUBRAND MD1C oil-free diaphragm vacuum pump. Compact and whisper-quiet, this pump is suitable for most benchtop rotary evaporators, and powerful enough for excellent results with high boiling point solvents.

Ultimate Vacuum: 1.5Torr, (2mbar)
Free Air Capacity: 25lpm, (0.88cfm)

The MD1C oil-free vacuum pump is significantly quieter than rotary vane pumps, waterjet aspirators, and competitive diaphragm pumps. Operation produces decibel levels comparable to a quiet, face-to-face conversation. With flowrates as much as 40% higher at working vacuum level than competitive pumps with similar free air capacity ratings. Flowpaths constructed of materials with excellent corrosion resistance, such

Vacuum Pumps: Vacuubrand MD1C Vacuum Pump

as PTFE heads and diaphragms. It is also equipped with a gas ballast to provide consistent pump performance with condensed vapors. The Vacuubrand chemistry-design diaphragm pumps eliminate the misting, handling, and disposal of contaminated oil associated with rotary vane pumps.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc., has been based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com