

Mettler Toledo AB204



Mettler Toledo AB204 Basic Level analytical balance, 210g x 0.1mg

Rating: Not Rated Yet

Price

Sales price \$1,225.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

!!!SOLD!!!

This Mettler Toledo AB204 Basic Level analytical balance has a range of 210g x 0.1mg, RS232 interface, full range taring, piece counting, percent weighing, & dynamic weighing.

Very clean

Please call us at (650) 871-5707 to check availability of this product!