

## Mettler Toledo AT200



Mettler Toledo AT200 Professional Level analytical balance, 205g x 0.1 mg

Rating: Not Rated Yet

**Price**

Sales price \$2,450.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

### Description

This Mettler Toledo AT200 Professional Level analytical balance, 205g capacity x 0.1 mg readability has Fully Automatic Calibration Technology (FACT), meets all GLP and ISO 9000 requirements, has a built-in bidirectional interface, and has an automatic draft shield that can be opened or closed at the press of a button.

Equipped with a 3.125" x 3.125" square weighing pan and full range taring.

**SOLD!!**

