

## Mettler Toledo AX205 Analytical Balance



Mettler Toledo AX205 Professional Level analytical balance, 220g x 0.01mg

Rating: Not Rated Yet

**Price**

Sales price \$3,850.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

### Description

This Mettler Toledo AX205 Professional Level analytical balance has a range of 220g x 0.01mg, comes w/fully automatic calibration technology (proFACT) & configurable draft shield w/automated door, and backlit graphic display w/touch screen operation.

### **SPECIFICATIONS - AX205**

Maximum load: 220 g

Readability: 0.01 mg

Repeatability at full load: 0.03 mg

Repeatability at gross load: 0.015 mg

Linearity: 0.1 mg

Eccentric load deviation at 1/2 maximum capacity: 0.25 mg

Sensitivity offset:  $1.5 \times 10^{-6}$

Sensitivity temperature drift: 0.001%/C

Sensitivity stability: 0.001%/a

Weighing time typical: 12s

Interface update rate: 7 /s

Balance dimensions: 241X505X293

Usable height of draft shield: 240 mm

Size of weighing pan: 80 x 80

This unit is very clean & was well maintained by the previous owner.