


Mettler Toledo Excellence Plus XP 8002 S



Fully automatic internal calibration (proFACT®) comes with history function, initiated by time and/or ambient temperature change by 1.5°C (2.7°F). Model 11277-128 features DeltaRange® capabilities, which deliver 10 times the readability anywhere within th

Rating: Not Rated Yet

Price

Sales price \$2,250.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

!!!! - RECENTLY SOLD - !!!!

Operating menu can be viewed in English, German, French, Spanish, Italian, Russian, and Japanese. All models can download software, display various units of weight, tare over entire weighing range, display statistics with \pm and % display, count with reference optimization, perform formula weighing with security protection, calculate with factors, manage users with password protection, store a history of setting changes, schedule balance checks, and monitor for allowable minimum weight.

All models feature a stainless steel platform, a color display with touch screen operation, sensors for hands-free operation, and feedthrough for below-the-balance weighing. Peripherals can be connected via the built-in RS-232 interface or an optional second interface (purchased separately). The following interfaces are available: Bluetooth® (wireless), Ethernet, PS/2, USB, LocalCAN™, or MiniMettler. For operation from 5 to 40°C (41 to 104°F) and between 25% and 80% relative humidity. Electrical: 110/230V, 50/60Hz, 0.8A.

SOLD!!!!