

Mettler Toledo FiveGo pH Meter



Model F2, meter only, AC/DC input 230 V AC, universal plug set

Rating: Not Rated Yet

Price

Sales price \$175.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

Ideally suited for mobile applications in the laboratory, at-line, or outdoors, the new FiveGo portable instruments provide high quality pH/mV, conductivity or dissolved oxygen measurements with the simple click of a button. The robust and waterproof design makes the FiveGo instruments the perfect companion to get reliable measurement data even in harsh conditions, all at a reasonable price.

Properties

BNC/cinch connector

display 3.1 in. LCD

model F2

pH measuring range: 0-14 pH unit, accuracy: ± 0.01 , resolution: 0.01

temperature measuring range: 0-100 °C, accuracy: ± 0.5 °C, resolution: 0.1 °C

battery power: 4 AAA

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com