

Mettler Toledo FP62



Mettler-Toledo FP62 Melting Point Apparatus in perfect working condition. Barely used. Equipped with RS232 data interface. Determines melting point of substances in the range from room temperature to 300 C.

Rating: Not Rated Yet

Price

Sales price \$2,850.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

This Mettler-Toledo FP62 Melting Point Apparatus is in perfect working condition & is barely used. It determines melting point of substances in the range from room temperature to 300 C.

Specifications:

Temperature range: from room temperature, +3 C to 300 C.

Heating / cooling rate: from 0.1 C/min to 10 C/min. Readability: 0.1 C Reproducibility at heating rate, 0.2 C/min: ± 0.1 C

Absolute accuracy of melting point determination at 0.2 C/min:

Up to 200 C, ± 0.5 C

From 200 C to 300 C, ± 0.8 C

The required sample quantity is about 1 – 3 mg.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com