

Mettler Toledo Halogen Moisture Balance HR73



Temperature range: 50-200C; Maximum sample weight: 71g

Rating: Not Rated Yet

Price

Sales price \$2,500.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

Mettler Toledo Halogen Moisture Balance HR73

Temperature range: 50-200C

Maximum sample weight: 71g

Readability of the balance: 1mg

Readability of the result: .01%

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.

226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com