

## Mettler Toledo MS3002TS Top-loading Balance



Capacity x Readability: 3,200 g x 0.01 g

Rating: Not Rated Yet

**Price**

Sales price \$1,475.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

### Description

The Mettler Toledo MS3002TS precision balance has a capacity of 3,200 grams with a readability of 0.01 grams and a pan size of 6.7 x 7.9 inches (170 x 200 mm).

### Technical Specifications:

Capacity: 3,200 g

Readability: 0.01 g

Pan Size: 6.7 x 7.9" (170 x 200 mm)

## Balances: Mettler Toledo MS3002TS Top-loading Balance

---

Repeatability: 0.01 g

Linearity deviation: 0.02 g

Typical Values

Repeatability: 0.007 g

Linearity deviation: 0.006 g

Sensitivity offset (test weight): 0.018 g (3,000 g)

USP minimum sample weight(5% load, k=2, u=0.1%): 14 g

Minimum sample weight (5% load, k=2, u=1%): 1.4 g

Settling time: 1.5 second

Interface: sUSB (device), USB (host), RS232, Bluetooth (optional)

Embedded applications: Weighing, Statistics for all applications, check weighing, totaling,

piece counting, formulation, percent weighing, factor weighing,  
dynamic weighing, density, dosing, PC Direct (for easy data transfer)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)