


## Mettler Toledo PB3002-S DeltaRange Balance



Readability: 0-600 grams: 0,01g; 600-3100 grams: 0,1g

Rating: Not Rated Yet

**Price**

Sales price \$925.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

**Description**

Balance from Mettler, type PB3002-S DeltaRange.

Max 3100 g; Min .5 g

e= 0.1 g; d= 0.01 g / 0.1 g

**Technical specifications:**

Voltage: 115V, 50Hz

Power: 6W

Readability: 0-600 grams: 0,01g; 600-3100 grams: 0,1g

Dimensions weighing platform (Ø): 18 cm

External dimensions (W) x (D) x (H): 24 x 32 x 9 cm

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)