


Mettler Toledo PB3002-S Precision Balance



Readability: 0.01 g; Max. capacity: 3100 g

Rating: Not Rated Yet

Price

Sales price \$850.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

Precision balance.

Specifications:

Readability: 0.01 g

Max. capacity: 3100 g

Repeatability (s): 0.01 g

Linearity: $-/+ 0.02$ g

Sensitivity drift: 6 ppm/ °C

Typical stabilization time: 2.0 s

Balances: Mettler Toledo PB3002-S Precision Balance

Adjustment weight: 2000 g

External dimensions of balance (W/D/H): 245/321/89 mm

Weighing pan \varnothing : 180 mm

Max. height above weighing pan: 141 mm

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com