

Mettler Toledo PB602-S



Mettler Toledo PB602-S toploading balance, 610g x 0.01g

Rating: Not Rated Yet

Price

Sales price \$675.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

The Mettler Toledo PB602-S - toploading balance is microprocessor controlled, large 7" diameter stainless steel weighing platform, built-in RS-232 interface, piece counting, percent weighing, dynamic weighing, and full range taring. 610g capacity by 0.01

repeatability. Very clean.