

## Mettler Toledo RB04-50DS MultiMax Automatic Reactor



System includes RB04-50DS Reactor Box, Docking Station, Power Box, Electronic Box, Universal Control Box, two Dispenser Boxes with 1000 $\mu$ L and 10mL syringes, Stirrer Control Set, four reactors of 50mL volume, four 160-1500 rpm stirrer motors, four InLab102

Rating: Not Rated Yet

### Price

Sales price \$28,000.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

### Description

System includes RB04-50DS Reactor Box, Docking Station, Power Box, Electronic Box, Universal Control Box, two Dispenser Boxes with 1000 $\mu$ L and 10mL syringes, Stirrer Control Set, four reactors of 50mL volume, four 160-1500 rpm stirrer motors, four InLab1022 pH probes, four Tr sensors, USB-4COM-SI-M RS-422/485 USB to serial adapter, also comes with a Dell Latitude D810 laptop computer operating on Windows XP OS with Mettler WinRC MM v2.12 Reactor Control software installed and with the Extended Calorimetry option, 115V, 60Hz, 250VA.

[Biotech Equipment Sales, Inc.](#)  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)