

Mettler Toledo Top Loading Balance SB12001



Maximum Capacity: 12,100,000mg, 12100g

Rating: Not Rated Yet

Price

Sales price \$975.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

Description

Sturdy lab balance. All-metal digital lab scale with capacity of up to 12100g (12,100,000mg) with a readability of 0.1g (100mg). Toggles between 14 weight units, and has useful functions such as piece counting, percent weighing, and dynamic weighing with either automatic or manual start. Settings are easily accessed from the up-front digital control panel with a backlit display that is easy on the eyes. Removable 9 x13 inch pan offers easy cleaning. With a low-profile design, it's easy to load and tare. Includes Monoblock weighing technology inside and an RS232 interface used to hook up a printer if needed.

Specifications:

Maximum Capacity: 12,100,000mg, 12100g

Readability: 100mg, 0.1 g

Linearity: ± 300 mg, 0.3 g

Power: 8-14.5V 50/60Hz 5Va, 9.5-20V, 5W

Dimensions: 13 x 15 x 4"

Pan Dimensions: 9 x 13"

Weight: 16lbs

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com