

## Mettler Toledo XPR404S Analytical Balance



Capacity of 410 g; readability of 0.0001 g; Ø2.7 (68 mm) pan size.

Rating: Not Rated Yet

**Price**

Sales price \$5,350.00

[Ask a question about this product](#)

Manufacturer [Mettler Toledo](#)

### Description

Mettler-Toledo XPR small platform analytical balances deliver outstanding weighing performance and support the highest requirements for data integrity and compliance.

**Technical Specifications:**

Maximum Capacity: 410 g

Readability: 0.1 mg

Repeatability (typical): 0.06 mg

Minimum Weight (USP), Typical: 120 mg

Weighing Pan Dimensions (DxW): 68 x 68 mm (SmartPan)

## Balances: Mettler Toledo XPR404S Analytical Balance

---

Linearity:  $\pm 0.2$  mg  
Repeatability (Test Weight): 0.1 mg (400 g)  
Settling Time: 2 seconds  
Dimensions (DxHxW): 411 x 368 x 214 mm  
Adjustment: Internal / proFACT Advanced  
Linearity: 0.2 mg  
Material Number(s): 30315628

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**

**226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**