

## Microfluidizer LM20, 2024 model -- only used 10 times!



Up to 30,000 psi (2,068 bar)

Rating: Not Rated Yet

**Price**

Sales price \$23,500.00

[Ask a question about this product](#)

### Description

The Microfluidizer LM20 (2024 model, manufactured by Microfluidics / Quadro Engineering) is a compact, electrically powered, high-shear fluid

processor / high-pressure homogenizer designed for laboratory-scale processing of small sample volumes. It uses intense shear forces, impact, and cavitation in a fixed-geometry chamber to create uniform submicron emulsions, dispersions, liposomes, and achieve efficient cell disruption with consistent, repeatable results. It's a premium lab "workhorse" for nanotechnology-level processing.

#### Key Features

**High Operating Pressure** — Up to 30,000 psi (2,068 bar), enabling superior particle size reduction and cell lysis at lower pressures compared to many alternatives.

**Small Sample Capability** — Minimum sample size of ~14 ml (with low hold-up volume), making it perfect for precious or limited materials.

**Digital Control & Consistency** — Siemens touchscreen HMI with digital pressure control, stroke counter, and maintenance reminders for precise, repeatable processing.

**Efficient Cooling & Temperature Control** — Designed to limit sample temperature rise during processing (max feed temp 165°F / 75°C), helping preserve heat-sensitive biomolecules.

**Scalable & Robust Design** — Benchtop unit (approx. 23" H x 33" W x 26" L, ~230–250 lbs) with linear scalability to larger Microfluidizer models; built with high-quality components including a 300 ml glass reservoir and diamond or ceramic Interaction Chamber.

Manufactured in Canada by Quadro Engineering

#### Other Specifications

Model: LM20

1100 PSI accumulator charge

Power: 208-230/60Hz or 190-220/50Hz, 3-phase

Weight: 230 lbs (104 kg)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)