

Multitron Floor Model Incubator Shaker Triple Stack with CO2/Humidity



Manufactured in 2021.

Rating: Not Rated Yet

Price

Sales price \$72,500.00

[Ask a question about this product](#)

Description

Manufactured in 2021, Custom configuration including x2 Multitron chambers (top/middle) configured for heating, shaking, CO2 control, humidity control, and both with a 25mm (1in.) orbital throw, bottom chamber is configured for heating, shaking, CO2 control, humidity control and has a 3mm orbital throw intended for use with microplates, 120V.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com