

## New Brunswick C-25KC



Refrigerated Incubator Shaker, large capacity floor model with refrigeration unit that provides temperature control from 5°C for below ambient investigations to 60°C

Rating: Not Rated Yet

**Price**

Sales price \$4,950.00

[Ask a question about this product](#)

Manufacturer [New Brunswick Scientific](#)

### Description

**!!! SOLD !!!**

This New Brunswick C-25KC Refrigerated Incubator Shaker is a large capacity floor model with refrigeration unit that provides temperature control from 5°C for below ambient investigations to 60°C. It's hermetically sealed, w/a large capacity compressor coupled with high velocity air circulation provides fast cooling response and precise temperature. Rotates from 40 to 400 rpm. Comes w/a universal platform & was manufactured in 2004 w/only 5900 hours of run time.

Very nice!

**Please call us at (650) 871-5707 to check availability of this product!**