

New Brunswick Innova 2100 Platform Shaker



New Brunswick Scientific Innova 2100 Platform Shaker

Rating: Not Rated Yet

Price

Sales price \$1,800.00

[Ask a question about this product](#)

Manufacturer [New Brunswick Scientific](#)

Description

Operating Features

- Maximum flask size: 6000mL
- Shaking speed: 25 to 500rpm
- Shaking orbit: Circular, 0.75 in. dia. (19mm)
- Tactile keypad with digital LED display for speed, time remaining, and elapsed time
- 0.1 to 99.9-hour timer
- Speed and timer alarms

Design Features

- Triple-action Gyrotory drive delivers smooth, uniform shaking
- Microprocessor stabilizes rate to 1rpm of setpoint
- Feedback loop compensates for fluctuations in voltage, load imbalance or workload
- Low-profile design; can be stacked inside environmental chamber
- Retains setpoints in memory; automatic restart if power is interrupted
- Acceleration circuit prevents sudden starts and minimizes splashing

Electrical requirement: 100/120/230/240V, 50/60Hz, 80VA, single phase

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office

650-276-7487 fax

sales@biotechequipmentsales.com