

## New Brunswick Innova 4230



Tall form New Brunswick Innova 4230 Incubator Shaker offers refrigeration as well, temperature range 20 degrees below ambient- minimum setpoint of 4C, to 80C, with shaking 18 x 18 platform & shelves above for standard incubation

Rating: Not Rated Yet

**Price**

Sales price \$4,250.00

[Ask a question about this product](#)

Manufacturer [New Brunswick Scientific](#)

Description

**!!! SOLD !!!**

This tall form New Brunswick Innova 4230 Incubator Shaker offers refrigeration as well, temperature range 20 degrees below ambient- minimum setpoint of 4C, to 80C, with shaking 18 x 18 platform & shelves above for standard incubation.

**Please call us at (650) 871-5707 to check availability of this product!**