

New Brunswick Scientific C24 Incubator/Shaker



The New Brunswick Scientific C24 is a benchtop incubator shaker with a temperature range of Ambient +7 to 60 deg C. 120V 50/60 hz. 5 x 2 Capacity. Utility Tray Platform.

Rating: Not Rated Yet

Price

Sales price \$2,650.00

[Ask a question about this product](#)

Manufacturer [New Brunswick Scientific](#)

Description

SOLD!!!

- **Temperature Range:**
E24: 7°C above ambient to 60°C.
E24R: 15°C below ambient to 60°C.
Minimum setpoint 4°C
- **Agitation:** 50 - 400 rpm, ± 2 rpm
- **Orbit:** 3/4" (1.9 cm) diameter orbit
- **Uniform Temperature Control** provided by self-correcting microprocessor feedback control and tightly constructed chamber. High-velocity fans ensure rapid thermal recovery
- **Temperature & Speed Can Be Calibrated** via keypad
- **Crystal Clear Lid Provides Excellent Visibility** minimizing the need to open the chamber; lifts out of the way for easy accessibility to samples
- **RS-232 Provided for Data Logging**
- **Safety Features** include automatic cut out of shaker when lid is opened to protect operator; and fail-safe thermostat that turns off heat if high-temperature limit is exceeded to protect cultures
- **Includes Utility Tray Platform**
- **Electrical Requirement** - 100/120V, 50/60Hz, 945VAC

Please call us at (650) 871-5707 to check availability of this product!