

Nuair Biosafety Cabinet Class II 425-600

PRODUCT
Image
Unavailable

PRODUCT
Image
Unavailable

The Nuair 425 LabGard Class II A2 biological safety cabinets give you the option of connecting it to a facility HVAC system, or exhausting it back into the room with internal HEPA Filters.

Rating: Not Rated Yet

Price

Sales price \$4,250.00

[Ask a question about this product](#)

Manufacturer [Nuair](#)

Description

The cabinets airflow is 30% exhausted / 70% recirculated to minimize cross-contamination of low to moderate risk biologicals in the absence of volatile toxic chemicals. The cabinet control panel includes a motor voltage regulator, audible window alarm, differential pressure gauge, power switches for exterior mounted fluorescent and ultraviolet lights. Also, interior outlets and a blower motor switch are also included.