

## Panasonic MCO-18ACL-PA CytoGrow Doublestack CO2 Incubator



Temperature range: +5°C above ambient to +50°C

Rating: Not Rated Yet

**Price**

Sales price \$6,850.00

[Ask a question about this product](#)

Manufacturer [Panasonic](#)

### Description

The Panasonic MCO-18ACL-PA CytoGrow CO2 Incubator was designed with passive contamination control.

### Specifications:

Temperature range: +5°C above ambient to +50°C

Outer door: painted steel

Access port: 1: back side (1.2" [30 mm] diameter)

## Incubators: Panasonic MCO-18ACL-PA CytoGrow Doublestack CO2 Incubator

---

Insulation: rigid polyurethane foamed-in-place  
Temperature control range: ambient temperature +5°C to +50°C  
Temperature uniformity:  $\pm 0.25^\circ\text{C}$  (ambient temperature: 25°C, setting: 37°C, CO2: 5 %, no load)  
Max. power consumption: Max 300 W  
CO2 control range: 0% to 20%  
CO2 controller: Microcomputer, Input by keypad  
CO2 gas connection: 0.157 inch to 0.36 inch (4 mm to 6 mm) diameter tube can be connected  
CO2 gas pressure: 0.03 MPa(G) (0.3kgf/cm2(G), 4.3 psi(G))  
CO2 variation:  $\pm 0.15\%$  (ambient temperature: 25°C, setting: 37°C, CO2: 5%; no load)  
CO2 sensor type: Thermal Conductivity  
Contamination controls: InCu-SaFe copper enriched stainless steel interior and shelves, optional SafeCell UV light  
Display: digital softkeys  
Heating system: direct hear and air jacket system  
Humidifying system: natural evaporation with humidifying pan, 95% R.H.  $\pm 5\%$  R.H.  
Interior surfaces: InCu-saFe stainless steel with copper enriched alloy  
Shelves: 17.7" x 17.7" (450 x 540 mm)  
Effective capacity: 6.0 cu. ft. (170 L)  
Interior volume: 6.0 cu. ft. (170 L)  
Interior Dimensions (W x D x H): 19.3" 20.6" x 26.2" (490 x 523 x 665 mm)  
Weight: 203 lb. (93 kg)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)