

## Perkin Elmer FlexDrop III EX Precision Reagent Dispenser



Designed for rapid dispensing of reagents, buffers, solvents and cells for high throughput assay preparation and compound storage.

Rating: Not Rated Yet

**Price**

Sales price \$11,000.00

[Ask a question about this product](#)

Manufacturer [Perkin Elmer](#)

### Description

1-2 reagent capacity, 200nL-2mL dispense volume range in 0.1 $\mu$ L increments, compatible with SBS standard 96-, 384- and 1536-well microplates in shallow and deep well formats, 2x8 precision valves dispense head configuration, offers extra capacity using an integrated 50-plate stacking mechanism, user friendly programming via membrane keypad for up to 99 dispense methods, integration capabilities with other devices and workstation platforms via RS-232 interface, comes with two 50-plate stackers, 100-130V, 50/60Hz, 60VA max.

[Biotech Equipment Sales, Inc.](#)  
[226 Miller Ave.](#)  
[South San Francisco, CA 94080](#)  
[650-871-5707 office](#)  
[650-276-7487 fax](#)  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)