

PHCbi Cell-IQ MCO-170AICUVHL-PA Dual-Chamber Air-Jacketed CO2 Incubators



Comes with 3 shelves per chamber and rolling dolly cart

Rating: Not Rated Yet

Price

Sales price \$8,250.00

[Ask a question about this product](#)

Description

5.83 cu.ft. (165-Liter) volume per chamber, +5°C above ambient to 50°C temperature range, features direct heat and air jacket system, PID control for CO2 and temperature, dual IR CO2 sensors, Incusafe copper enriched stainless steel interior, Safe Cell UV light, touch screen panel, H2O2 vapor sterilization / decontamination function, comes with 3 shelves per chamber and rolling dolly cart, 110-120V, 60Hz.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com