

## PolyScienceWBE28 General Purpose Water Bath, 28L



Displays actual and set point temperature simultaneously

Rating: Not Rated Yet

**Price**

Sales price \$675.00

[Ask a question about this product](#)

Manufacturer [Polyscience](#)

### Description

General Purpose Water Bath 28 Liter

Model #: WBE28

#### Key Features:

Displays actual and set point temperature simultaneously

Programmable pre-sets for frequently used temperatures

Integrated timer

Hinged gable cover

Calibration off-set feature

Reservoir size is 28 liters

Reservoir made of 304 Stainless Steel

Reservoir drain

Primary and automatic safety thermostats

Also compatible with waterless bath beads

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)