


## Qiagen TissueLyserII



with 24 place adaptors.

Rating: Not Rated Yet

**Price**

Sales price \$6,500.00

[Ask a question about this product](#)

### Description

Research grade sample preparation system

An efficient solution for high-throughput sample preparation. This system quickly and effectively disrupts and homogenizes biological samples, with an emphasis on preventing cross-contamination. It's versatile, capable of processing diverse sample types, and integrates well with QIAGEN's tissue management workflow. The TissueLyser II releases high-quality DNA, RNA, and proteins, aiding further purification and analysis. It's an essential tool for enhancing biomolecule yields in genotyping, gene expression, and proteomics applications.

### Specifications:

Capacity: up to 192 samples simultaneously, depending on the adaptors used

Sample Size: 1.2ml collection microtubes

Oscillation Frequency: 3 - 30 Hz (180 - 1800 oscillations per minute)

Sample Disruption Type: bead mill

Power Requirements: 120V 15A

Voltage: 100-240V

Frequency: 50/60Hz

Power: 150W

Fuse: T; 3.15A (2x)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)