

## Rainin EDP 3 200ul 12 channel Electronic Pipettor



Includes Li-ion battery.

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [Rainin](#)

### Description

Includes Li-ion battery. Increased productivity. Ultra-light one-hand tip ejection. Autopace for auto-dispense plus versatile modes. Uses LTS tips. Order charging devices separately.

Manufacturer number:	E3-200
Brand:	EDP 3 Electronic Pipette
Manufacturer:	Rainin
Type:	Single Channel Electronic Pipette
Volume:	20-200 ul
Display:	Electronic
Ejector:	Yes
Autoclavable :	No
Speed Control:	Yes
Programmable:	Yes

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**