

Rainin Pipet-Lite L5000 Magnetic-Assist Pipettor

PRODUCT
Image
Unavailable

PRODUCT
Image
Unavailable

Pipet-Lite: Reduce total force by 70%. Discover hand-friendly pipetting. Magnetic assist helps locate and hold zero position for improved results. LTS LiteTouch Tip Ejection for ultralight tip ejection. Reduce your risk of RSIs. SL models uses standard un

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Rainin](#)

Description

Pipet-Lite: Reduce total force by 70%. Discover hand-friendly pipetting. Magnetic assist helps locate and hold zero position for improved results. LTS LiteTouch Tip Ejection for ultralight tip ejection. Reduce your risk of RSIs. SL models uses standard universal tips, L models uses LTS tips only.

Part number:	L-5000R
Manufacturer number:	L-5000
Brand:	Pipet-Lite Magnetic-Assist Pipette
Manufacturer:	Rainin
Type:	Single Channel Manual Pipette
Volume:	500-5000 ul
Display:	Digital
Ejector:	Yes
Autoclavable :	No
Speed Control:	No
Programmable:	No