

## REVCO ULT2586-9-A31 Ultima II -80 freezer



Dimensions: 77.5 inches high by 48 inches wide by 36 inches deep. 25 cubic foot capacity.

Rating: Not Rated Yet

**Price**

Sales price \$5,500.00

[Ask a question about this product](#)

Manufacturer [REVCO](#)

### Description

Revco Ultima II -86 Freezer.

Designed for clinical materials storage, these sophisticated freezers provide high security for laboratory samples. Units feature the InterLogic microprocessor control system, which is simple to use, yet sophisticated in performance with direct digital setpoint and display, keyed power and alarm switch, and touchpad data entry. The lifeguard compressor protection system monitors compressor performance and automatically adjusts operation. Automatic incident Monitor initiates a fast flash warning if the system has encountered problems and records and displays the minimum and maximum temperatures achieved. An extreme ambient alert provides an early warning when ambient temperature is high enough to affect compressor operating conditions. Units also feature voltage boost, surge suppressor, clean filter indicators, remote alarm controls, and RS-232 data port.

- Model ULT2586-9-A36
- 25 cu ft.
- Unit has 5 interior compartments each with an individual door.
- Temp range -50 to -86 deg C.
- Digital control for temp and alarms